



LERNAPHARM

Specialized in the Development and the Production of Antiseptics and Disinfectants

● **Over A Quarter Century Of Service**

● **International Presence**

**PHARMACEUTICAL
MANUFACTURING**



**Promoting
Health And Wellness**





LORIS™ Antiseptic Isopropyl Alcohol Pad

Product Code : 103-03

- Topical skin cleanser prior to injection
- Each pad saturated with 70% Isopropyl alcohol
- Individually wrapped
- Single-used
- Lint-free
- Latex-free
- External use only

Characteristics

- Effective in destroying harmful bacteria
- Helps to prevent the spread of bacteria
- Approved by Health Canada as an antiseptic skin cleanser
- Recommended by CDC for the skin preparation before central & peripheral venous catheter insertion, peripheral arterial catheter insertion, and during dressing changes



How To Use

Tear open the packaging. Rub skin with pad briskly in a circular motion from injection site outward. Allow to dry.



2323 Halpern
St-Laurent (Montreal)
Quebec, Canada, H4S 1S3

Tel: + 514-331-4634
Toll-Free: 1-866-322-4634








LORIS™ Antiseptic Isopropyl Alcohol Pad

Product Code : 103-03

Minimal Acceptance Criteria for Antiseptics Microbial Limit Test

Test	USP Nonsterile Dosage Forms (Cutaneous Use)	Lernapharm finished product
Total Aerobic Microbial Count (TAMC)	≤100 cfu / mL	<10 cfu / mL
Total Combined Yeast and Molds Count (TYMC)	≤10 cfu / mL	<2 cfu / mL
Staphylococcus aureus detection	Absent	Absent
Pseudomonas aeruginosa detection	Absent	Absent
B. Cepacia	Not specified	Absent

Packaging Information

Product	Quantity	Size (L,W,H)	Weight
Wipe 	-	65 mm x 30 mm	
Pouch 	1 x Wipe per Pouch	50 mm x 50 mm	
Display Box 	200 x Pouch per Display Box	110 mm x 108 mm x 56 mm	
Case 	20 x Display Box per Case	34 cm x 22 cm x 25 cm	4 kg
Pallet 	144 x Case per Pallet	122 cm x 102 cm x 201 cm	598 kg